



# YSI 2900C Biochemistry Analyzer

FDA CLASS II REGISTERED DEVICE



FDA 510(k) Registered Device

**YSI 2900C Biochemistry Analyzer is the flexible, reliable platform that supports 21CFR Part 11 compliance, and the accuracy and precision your work requires.**

Here's how 50 years of experience with patented enzyme electrode technology can support you:



Results in 60 seconds or less



Manual, single-sample loading,  
or full automation of up to 96 samples



Programmable, automated  
calibration



Unaffected by many of the physical and  
chemical interferences that challenge  
other methods



a xylem brand

# YSI 2900C Biochemistry Analyzer

Feature	Specifications
Certifications	ETL, RoHS
Compliance	FDA registered 510(k) - k210933
Connectivity / Communications	USB
Graphic Display	Color Touchpanel User Interface
Measurement Range	Glucose: 5-900 mg/dL
Memory	SSD Solid State Disk, data downloadable to flash drive
Operating Temperature	Operating Temperature Conditions: 15 - 35 °C and 10 - 75 % Relative Humidity
Parameters Measured	D-Glucose (Dextrose) in human whole blood, capillary fingerstick blood, plasma, serum, and aqueous solution
Power	100 - 240 VAC, 50 - 60 Hz
Precision	2% or 2.0 mg/dL, whichever is greater
Unit of Measure	mg/dL, mg/L, g/L, %, mmol/L
Warranty	1 year

SPECIFICATIONS XA00274



Xylem Lab Solutions  
 1725 Brannum Lane  
 Yellow Springs, OH 45387

+1.937.767.7241

info@ysi.com

ysi.com

© 2024 Xylem Inc. XA00274 1024



[YSI.com/2900C](https://www.ysi.com/2900C)